



Angela Safety

Angela minimizes the formation of free radicals
by remarkable reducing the protein which causes hypersensitivity.

Angela Protein Load : 11~12 μg/ml

 Angela has an excellent safety profile due to its very low levels of endotoxin.

Angela Endotoxin Unit Level: Max.0.25 Eu/g





Angela...

- Angela is a derivative of hyaluronic acid which is biocompatible with human body.
- Angela has an excellent physical property due to highly developed manufacturing technology.
- Angela is a biograft product which has a high complex viscosity & excellent resistance against breakdown enzymes.
- Angela has a brilliant safety profile for the body.
- **Angela** got advantages from both of bi-phasic Restylane(Restylane)'s and mono-phasic(Juvederm)'s.

Therefore, It has smooth soothing power, excellent volume, cohesiveness, duration, stability after injection.

Ex) Juvederm's softness & Restylene's volume

- Angela doesn't affect to swell after treatment (about 1 out of 10,000).
- Angela duration: about 12 months.



Premium HA

- Premium HA used to produce Angela, which certified by USA(FDA) and Europe(EDQM).
- HA from streptococcus zooepidemicus cultivation.

ANGELA FILLER LINE			
Product name	Sub - Q	Deep	Fine
Туре	For fine wrinkle	For deep wrinkle	For volume
Injection part	Side of eyes, ETC.	Smile lines, forehead, cheek, ETC.	Nose, chin, ETC.
Volume	1ml	1ml	1ml
HA concentration	24mg	24mg	24mg
Needle size	30G	28G	27G
Injection force	8.12	8.86	15.55





How to use

■ Guide line of parts to inject.



Sub-Q



Deep



Fine





Before / After

⊳Nose





•••

Before / After

⊳ Frown line



▷ Nasolabial Folds





Angela

Graft / Prosthesis / Biomaterial

